Mattress Bariatric

Anti-decubitus mattress with extra high comfort especially designed for use on bariatric beds. The mattress consists of 10 cm visco-elastic foam and 10 cm fire retardant comfort foam.

Core

- The core consists of:
 - An upper layer of 10 cm of high quality visco-elastic memory foam with unique thermo active properties and a density of 80-85 kg/m³. Moreover, the visco-elastic foam has an exceptional hardness of 3,5 kPa so the mattress volume does not decrease after heavy loads.
 - A highly elastic lower layer of 10 cm fire retardant polyurethane foam with a density of 45 kg/m³. The
 polyurethane foam is very elastic, fire retardant and thermosetting.
- Pressure relieving with optimum temperature sensitivity.
- Suitable for medium risk of decubitus.
- The mattress retains its elasticity even after extensive use.
- Bacteriostatic treatment against pathogenic and allergenic micro-organisms.
- Smooth and excellent flexibility for use on an adjustable bed.

Cover

- Bi-elastic cover (stretching chain: > 120%, warp: > 130%).
- The weight of the cover is 205 g/m².
- The material of the cover is completely adjusted to optimal pressure distribution and a low frictional resistance.
- The cover is waterproof, vapor permeable, breathable and urine resistant.
- The cover is antimicrobial, hypoallergenic and neutralizes pathogenic micro-organisms and allergens.
- Removable, with zipper on a long and short side for easy replacement.
- Washable at 95 °.



Dimensions

10604 | Bariatric mattress

- L 195 x B 90 x H 20 cm



Materials and finishing

- <u>Core</u>: Visco-elastic memory foam (80-85 kg/m³), polyurethane foam CMHR (45 kg/m³).
- Cover: Bi-elastic, polyurethane coating on a base of polyamide.
- The mattress and the cover meet fire standard CRIB 5.
- The mattress and cover are ecologically produced without CFCs. Neither the production nor the possible destruction processes are environmentally harmful.
- Withstands the common cleaning products.

For more information about the **materials**: please consult our material lists.

Options

- 6 filler mattresses (width: 10 cm)



- 6 filler mattresses (width: 20 cm)



