



# AVENTO

*high/low*

# Avento high/low

The Avento high/low relax chair is a real blessing for the patient. The chair is equipped with a high/low system and a relax mechanism that inclines the backrest, seat and legrest simultaneously into a "zero gravity" position. The patient enjoys an exceptional support in the Avento high/low chair.

Have a try,  
you'll love it!



## 01. HEIGHT ADJUSTMENT

*Comfort for everyone*

Height adjustment between 53 and 83 cm thanks to a hydraulic pump easy to operate by means of a pedal at both sides of the chair.

## 02. RELAX MECHANISM

*A blessing for the patient*

Less pressure points on the body and a good blood circulation thanks to the patented mechanism "Zero Gravity".

## 03. CENTRAL BRAKE

*Mobile and flexible*

The chair is extremely manoeuvrable and mobile, thanks to the large castors and the central brake system on 4 castors.

## 04. DESIGN

*Hidden functionality*

The mechanism is completely hidden by a design cover cap in synthetic material. The rounded corners of the cover cap make it more easy to clean the chair.

## 05. MODULARITY

*According to your needs*

The relax chair can be composed the way you want it. Various accessories make the chair more complete, such as a removable headrest, a dinner tray,...

## CONTACT US

*We are listening to you*

Haelvoet Belgium  
T +32 (0)51 48 66 95  
F +32 (0)51 48 73 19  
[info@haelvoet.be](mailto:info@haelvoet.be)

(\*) *Zero gravity position*: this seat/lying position is based on NASA data. A neutral body position converts all muscle power into a biomedical balance.  
(Source: "From outer space to you", Omni, March 1994)