
Accessories Hospital

The accessories we offer can be harmoniously connected with the existing cases and accessory holders. Haelvoet strongly pays attention to this, so that the accessories offer you the same top quality as our nursing beds. After all, only the best materials are good enough for you.



(*) Option





Colours

We would like to refer to our extensive colour gamut. Sample card + samples available on request.

(*) Option




Accessories care

The accessories we offer can be harmoniously connected with the existing cases and accessory holders. Haelvoet strongly pays attention to this, so that the accessories offer you the same top quality as our nursing beds. After all, only the best materials are good enough for you.

001 Lifting pole		006 IV holder	
	 Bent thick-walled steel tube. Handle with a manually adjustable ribbon. Automatic roll-up system (*). 		- Fixed IV holder, double-bent with 2 IV hooks Adjustable double-bent IV holder with 4 IV hooks (*).
002 IV holder for lifting pole		007 Flexible hand switch holder	
80707	- IV holder with 4 hooks.		Improves the accessibility of the hand switch for both the patient and nursing staff.
Dinner tray for sideguards		008 Raising aid	
	Made of a solid beech frame with a milled-in solid laminate tablet.Can be easily hung onto the foot panel when not used.		Optimal support when getting in and out of the bed. Ergonomic solid laminate design.
004 Accessory bracket		009 Horizontal urine bottle holder	
	 Chromium-plated bracket (310 mm) with 4 hooks. Installed in the accessory brackets (*). 		- Installed in the accessory hooks (*).
005 Bowl holder		010 Vertical urine bottle holder	
(4) Ootion	Installed in the accessory hooks (*). Possible to install 2 bowls.		- Installed in the accessory hooks (*).

(*) Option





Colours

We would like to refer to our extensive colour gamut. Sample card + samples available on request.

(*) Option

