

Advance Hospital

06293 The Advance is a very elegant, affordable and durable bedside table. The one-sided version is equipped with a refrigerator and an optional dinner tray.

One-sided version with refrigerator (storage cabinet with refrigerator, open compartment, drawer)

- 2 side panels, a bottom and a tabletop.
- The back is made of steel plate.
- The side panels are covered with 2 MDF posts with an aluminium look and a light grey melaminated chipboard of 19 mm thick.
- Noiseless absorption refrigerator of "Dometic Minicool - Silencio DS 400 BI".
- Capacity: 37 litres.
- Door with HPL finishing.
- Interior lighting with LED.
- Spiral cable.
- Metal design handle.
- An open compartment is accessible from on the front side.
- The drawer is made of steel plate and has a solid laminate front.
- The drawer glides on 2 synthetic rails with an end stop.
- The drawer is equipped with a removable synthetic drawer tray.
- The metal parts are protected by a light grey epoxy-coating.
- The tabletop has rounded solid laminate corners.
- 4 design castors Ø 65 mm.

Optional dinner tray Livia

- Folding.
- Easy to incline in both directions.
- The dinner tray has an upright edge.
- The folding part is fixed on a strong aluminium guiding profile that has been finished with an anodisation layer.
- The dinner tray is height-adjustable (between 79 and 112 cm) by means of a gas spring.
- Dimensions of the dinner tray: 63 x 40,5 cm.
- The gas spring is activated by means of a foot pedal.
- The dinner tray is made of fire-retardant ABS.
- The bedside table stands on 5 design castors Ø 65 mm.



Dimensions

06293 | with refrigerator | 49,5 x 49,5 x height 88 cm
with dinner tray Livia | 67 x 49,5 x height 88 cm



Materials and finishing

- Solid laminate
- Wood: MDF and melaminated chipboard (light grey).
- Metal: epoxy-coating (light grey).
- High-quality self-closing hinges.
- Withstands the common cleaning products.

For more information about the **materials, construction and maintenance**: please consult our material lists.

Options

- Livia dinner tray.
- Cylinder lock.
- Gallery.
- Cases for a telephone holder.
- Telephone holder.
- Removable tray with upright edges for the tabletop.
- 4 bumper castors.
- 4 castors with brakes.

Colours

We would like to refer to our extensive colour gamut.
Sample card + samples available on request.

