

Isolation bed

02465 The isolation bed is made of a solid welded metal frame in a wooden case. The bed is available with different mattress platform heights.


Wooden bed

- Solid welded metal frame made of tubes (35 x 35 x 2 mm).
- The wooden case around the metal frame is made of plywood panels with a lower indented plinth.
- The mattress platform is surrounded by an edge in solid beech with 5 attachment points for fixation straps at both sides. At the head and foot end, 1 attachment point for fixation straps is provided.
- The mattress platform is made of a removable plywood panel, thickness 23 mm.
- Rounded corners.



Dimensions

- Width: 90 cm. Length: 200 cm
- Height of the mattress platform: 35 / 45 / 65 cm
- Mattress platform dimensions: 200 x 90 cm
- Mattress dimensions: 195 x 85 cm
thickness: 12 to 14 cm



Materials and finishing

- Wood: Beech (solid), plywood.
- Metal: steel.
- Withstands the common cleaning products.

For more information about the **materials** and **construction**: please consult our material lists. We would like to refer to our extensive **colour gamut**.

Options

- Opening for emergency button.
- Non-removable mattress platform.
- Floor fixation.
- Mattress in fire-retardant foam with fire-retardant mattress cover.
- Adapted height of the mattress platform.

Colours

We would like to refer to our extensive colour gamut. Sample card + samples available on request.

