# Talo chair with castors

The Talo rest chair with castors is a comfortable and well-designed chair. The different foam density zones guarantee an excellent seat comfort.

# **Support**

- The support is made of 2 solid beech sides, connected to each other by means of a solid beech cross beam.
- The side parts are connected to each other by means of a glued mortise and tenon joint.
- The sides are firmly connected to each other by means of a solid beech cross beam and a metal frame.
- The ergonomic design of the armrests makes it easy to stand up or sit down.
- Comfortable rounded armrests of  $\pm$  7 cm.
- All edges and corners have been rounded.
- Four design twin castors Ø 100 mm with separate brakes.

### Seat

- Removable seat, made of a melaminated MDF board of 12 mm.
- The rounded seat in front stimulates a fluent blood circulation.
- Comfortable seat height of 47 cm.
- Upholstered with fire-retardant foam in 2 different density zones to guarantee an excellent seat comfort.
- The upholstery is finished with imitation leather and decorative stitching.

## **Backrest**

- Made of an ergonomically preformed plywood shell of 12 mm.
- Upholstered with fire-retardant foam in 3 different density zones to guarantee an excellent seat comfort.
- The upholstery is finished with imitation leather and decorative stitching.
- Large and ergonomic push bar made of synthetic material at the back.

We would like to refer to our extensive colour gamut. Sample card + samples available on request.

## **Dimensions**

14184 | W 68 x D 94 x H 124 cm

- Seat height: 47,5 cm
- Seat depth: 51 cm
- Seat width: 53,5 cm
- Height armrests: 65,5 cm





# Materials and finishing

- Wood: beech (solid), plywood, MDF.
- Metal: epoxy coating (light grey).
- Fire-retardant foam.
- <u>Upholstery</u>: see collection.
- Resistant to the common cleaning products.

For further information about used materials, constructions and maintenance: consult our materialfiles.

# **Options**

- Dinner tray attachment on armrests.
- Seat cushion with visco-elastic foam and a bielastic cover.
- Anatomical reduction cushion.
- Lateral headrest.
- Backrest with removable headrest.
- Closed side panels.
- Accessory brackets.
- Restraint bracket.
- Case for IV rod.
- Dinner tray support on backrest.



Colours