

Stool_{wood}

00953 Stable, simple and space-saving.
Ideal for visitors in a patient room.

Support

- The 4 legs are made of laminated wood. They are firmly screwed on the round seat Ø 32,5 cm.
- The wooden parts are finished with a coat of varnish.
- The wooden stool is very stackable.
- The legs are equipped with synthetic glides.


Seat

- The upholstered round seat is made of chipboard and has rounded edges.
- If it is not upholstered, then the seat is made of high-quality plywood and finished with veneer.
- The seat is upholstered with fire-retardant foam (thickness: 20 mm; density: 50 kg/m³) and finished with imitation leather (in case the customer chooses the upholstered version).



Dimension

00953 | 39,5 x 46 cm seat height



Materials and finishing

- Wood: beech (Laminated), plywood, chipboard.
- Metal: epoxy-coating (light gray).
- Fire-retardant foam.
- Upholstery: see collection.
- Resistant to the common cleaning products.

For further information about used **materials**, **constructions** and **maintenance**: consult our materialfiles.

Colours

We would like to refer to our extensive colour gamut.
Sample card + samples available on request.

